



SterileRight

Packaging Mfg. Inc.



Hospital & Dental

STERILIZATION PACKAGING



EN868/ISO11607/11140/13485/ GMP www.sterileright.com.tw




SterileRight


SterileRight

Packaging Mfg. Inc.

Your reliable
OEM & ODM
partner



● Flat Rolls



● Gusseted Rolls

● Self-seal Pouches



Features

- Medical grade paper 70gsm/ Multilayer non-tearing green film
- Medical grade paper 60gsm/ Multilayer non-tearing blue film
- Environmentally friendly printing / Multi-lined seals
- Higher performance of seal strength & clean peeling
- Printed with Steam & EtO indicators (Formaldehyde indicator upon request)
- Customized sizes and specifications supplied upon request

Sterilization methods

Suitable for **Steam**, **EtO** and Formaldehyde




See-through Sterilization Rolls & Pouches

Flat Rolls:

	ITEM NO.	SIZES	R/CTN	ITEM NO.	SIZES	R/CTN
	Medical grade paper / Multilayer non-tearing Blue film			Medical grade paper / Multilayer non-tearing Green film		
60 GSM	01.6001	50mmx200M	16	70 GSM	01.7001	50mmx200M
	01.6002	75mmx200M	8		01.7002	75mmx200M
	01.6003	100mmx200M	8		01.7003	100mmx200M
	01.6004	150mmx200M	4		01.7004	150mmx200M
	01.6005	200mmx200M	4		01.7005	200mmx200M
	01.6006	250mmx200M	2		01.7006	250mmx200M
	01.6007	300mmx200M	2		01.7007	300mmx200M
	01.6008	350mmx200M	2		01.7008	350mmx200M
	01.6009	400mmx200M	2		01.7009	400mmx200M
	Max. Size : 800mm x 200M			Max. Size : 800mm x 200M		

Gusseted Rolls:

	ITEM NO.	SIZES	R/CTN	ITEM NO.	SIZES	R/CTN
	Medical grade paper / Multilayer non-tearing Blue film			Medical grade paper / Multilayer non-tearing Green film		
60 GSM	02.6002	75mmx35mmx100M	8	70 GSM	02.7002	75mmx35mmx100M
	02.6003	100mmx40mmx100M	8		02.7003	100mmx40mmx100M
	02.6004	150mmx50mmx100M	4		02.7004	150mmx50mmx100M
	02.6005	200mmx55mmx100M	4		02.7005	200mmx55mmx100M
	02.6006	250mmx60mmx100M	4		02.7006	250mmx60mmx100M
	02.6007	300mmx70mmx100M	2		02.7007	300mmx70mmx100M
	02.6008	350mmx80mmx100M	2		02.7008	350mmx80mmx100M
	02.6009	400mmx80mmx100M	2		02.7009	400mmx80mmx100M
		Max. Size : 600mm x 100M			Max. Size : 600mm x 100M	



See-through Self-seal Pouches

Assured sterility / Protection seal / Thumb-notch
User-friendly pre-folded self-sealing



Accurate & Fast seal

	ITEM NO.	SIZES	BOXES/CTN	PCS/CTN	ITEM NO.	SIZES	BOXES/CTN	PCS/CTN
	Medical grade paper / Multilayer non-tearing Blue film				Medical grade paper / Multilayer non-tearing Green film			
60 GSM	03.6001	57mm x 133mm	90	18000	70 GSM	03.7001	57mm x 133mm	90
	03.6002	90mm x 162mm	30	6000		03.7002	90mm x 162mm	30
	03.6003	70mm x 257mm	40	8000		03.7003	70mm x 257mm	40
	03.6004	90mm x 257mm	30	6000		03.7004	90mm x 257mm	30
	03.6005	135mm x 283mm	15	3000		03.7005	135mm x 283mm	15
	03.6006	180mm x 335mm	15	3000		03.7006	180mm x 335mm	15
	03.6007	190mm x 358mm	10	2000		03.7007	190mm x 358mm	10
	03.6008	300mm x 380mm	5	1000		03.7008	300mm x 380mm	5
	03.6009	300mm x 474mm	4	800		03.7009	300mm x 474mm	4

Sterilization Tyvek® Rolls & Pouches

Features

- Made from uncoated Tyvek® 1073B, 1059B, 40L & high functional film
- Environmentally friendly printing
- Printed with Plasma indicator
- Standards: EN 868, ISO 11607, 11140
- Customized sizes and specifications supplied upon request

Sterilization methods

Suitable for Hydrogen Peroxide (Plasma), EtO, Gamma, and E-beam sterilization

1073B

1059B

Sterilization Tyvek® Rolls (1073B/Clear Film)

ITEM NO.	SIZES	R/CTN
01.7301	50mm x 100M	24
01.7302	75mm x 100M	16
01.7303	100mm x 100M	12
01.7304	150mm x 100M	8
01.7305	200mm x 100M	8
01.7306	250mm x 100M	4
01.7307	300mm x 100M	4
01.7308	350mm x 100M	4
01.7309	400mm x 100M	4

Max. Size : 800mm x 100M

Sterilization Tyvek® Rolls (1059B/Clear Film)

ITEM NO.	SIZES	R/CTN
01.5901	50mm x 100M	24
01.5902	75mm x 100M	16
01.5903	100mm x 100M	12
01.5904	150mm x 100M	8
01.5905	200mm x 100M	8
01.5906	250mm x 100M	4
01.5907	300mm x 100M	4
01.5908	350mm x 100M	4
01.5909	400mm x 100M	4

Max. Size : 800mm x 100M

Miscellaneous properties of Tyvek®

Microbial Barrier

Log Reduction Value(LRV)

5.2

Comparable Test Method : ASTM F1608

Bendtsen Air Permeability

mL / min

572

Comparable Test Method : ISO 5636-3

Spencer Puncture

J / m²

8756

Comparable Test Method : ASTM D3420

Microbial Barrier

Log Reduction Value(LRV)

4.7

Comparable Test Method : ASTM F1608

Bendtsen Air Permeability

mL / min

671

Comparable Test Method : ISO 5636-3

Spencer Puncture

J / m²

6830

Comparable Test Method : ASTM D3420



SterileRight

SterileRight Packaging Mfg. Inc.

Sterilization Tyvek® Pouches

Sterilization Tyvek® Rolls



40L

Sterilization Tyvek® Rolls (40L/Clear Film)

ITEM NO.	SIZES	R/CTN
01.4001	50mm x 100M	24
01.4002	75mm x 100M	16
01.4003	100mm x 100M	12
01.4004	150mm x 100M	8
01.4005	200mm x 100M	8
01.4006	250mm x 100M	4
01.4007	300mm x 100M	4
01.4008	350mm x 100M	4
01.4009	400mm x 100M	4

Max. Size : 800mm x 100M



2FS

Sterilization Tyvek® Rolls (2FS/Clear Film)

ITEM NO.	SIZES	R/CTN
01.5801	50mm x 100M	24
01.5802	75mm x 100M	16
01.5803	100mm x 100M	12
01.5804	150mm x 100M	8
01.5805	200mm x 100M	8
01.5806	250mm x 100M	4
01.5807	300mm x 100M	4
01.5808	350mm x 100M	4
01.5809	400mm x 100M	4

Max. Size : 800mm x 100M



Microbial Barrier

Log Reduction Value(LRV)

2.6
(Ref)

Comparable Test Method : ASTM F1608

Bendtsen Air Permeability

mL / min

2350

Comparable Test Method : ISO 5636-3

Spencer Puncture

J / m²

4553

Comparable Test Method : ASTM D3420

Microbial Barrier

Log Reduction Value(LRV)

3.2
(Ref)

Comparable Test Method : ASTM F1608

Bendtsen Air Permeability

mL / min

600

Comparable Test Method : ISO 5636-3

Spencer Puncture

J / m²

4903

Comparable Test Method : ASTM D3420

Sterilization Indicator Tapes

Distinguishes processed from unprocessed items
Adheres well to woven and non-woven materials

Item No. M901.18.0101

Steam Tape / Color : White or Green changed to Black
Size : 19mm x 50M / Roll

Item No. M901.18.0201

EtO Tape / Color : Brown changed to Green
Size : 19mm x 50M / Roll



Non-Woven Sterilization Wraps & Medical Crepe Paper

Available in a variety of sizes

Ideal package material during the Steam and EtO
sterilization cycle

Colors can be green, blue and white

Available width: 300mm x 300mm ~ 1200mm x 1200mm

Medical Crepe Paper



Non-Woven Sterilization Wraps



SterileRight

SterileRight

Packaging Mfg. Inc.

Fairness Efficiency Innovation Development

- **Design**
(Standards: EN 868, ISO 11607, ISO 11140)
- **Printing**
(Environmental friendly printing)
- **Pouch Making**
(Standards: EN 868, ISO 11607, ISO 11140)
- **Quality Inspection**
(Standards: EN 868, ISO 11607, ISO 11140)



Contact:

- 📞 **Office:** 886-2-2697-5499 **Fax:** 886-2-2697-5497
8F-10, No. 99, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City 221, Taiwan
- Factory:** 886-2-2455-5498 **Fax:** 886-2-2455-5497
1F, No. 33-7, Dahua 2nd RD., Qidu Dist., Keelung 206, Taiwan
- ✉ **E-mail:** service@sterileright.com.tw
- 🌐 **Website:** www.sterileright.com.tw