



# SterileRight

Packaging Mfg. Inc.



## **Medical Device** STERILIZATION PACKAGING



EN868/ISO11607/11140/13485/ GMP [www.sterileright.com.tw](http://www.sterileright.com.tw)



# Fairness Efficiency Innovation Development

- p04----- Reinforced Header Bags
- p05----- Peel Pouches
- p06----- FFS-Form Fill Seal Packaging
- p07----- Coated Tyvek® Lid
- p08----- Pre-made Pouches  
Breathable Strip Bags
- p09----- Aluminum Foil Bags  
Aluminum Foil / Tyvek® Header Bags
- p10----- Surgical Equipment Drape Packs

# Medical Device Packaging

## Reinforced Header Bags

### Features

Made from un-coated or coated Tyvek® and high functional film  
 Ideal for **heavy loading** surgery and procedure kits  
 Suitable for bulk sterile kits, drapes and gowns  
 Higher seal strength & superior breathability  
 Impact strength, puncture and tear resistance

Customized sizes and specifications supplied upon request

### Sterilization methods

Suitable for EtO and Irradiation



Burst Test

Patent in Taiwan



## Peel Pouches

### Features

Supplied with paper/paper, paper/film, Tyvek® /film and film/film  
 Chevron seal, thumb notch, with or without coating produced according to customer's need

### Sterilization methods

Suitable for EtO or Irradiation



Sterilization Tyvek® Pouch



Customized printing, sizes and specifications supplied upon request

Authorized converter of DuPont™ Tyvek®  
 Tyvek® is DuPont's trademark.

| Item No. | Sizes                           |
|----------|---------------------------------|
| 05.0101  | 387mm x 425mm (15.25" x 16.75") |
| 05.0102  | 413mm x 502mm (16.25" x 19.75") |
| 05.0103  | 438mm x 578mm (17.25" x 22.75") |
| 05.0104  | 464mm x 540mm (18.25" x 21.25") |
| 05.0105  | 476mm x 413mm (18.75" x 16.25") |
| 05.0106  | 514mm x 578mm (20.25" x 22.75") |
| 05.0107  | 514mm x 629mm (20.25" x 24.75") |
| 05.0108  | 546mm x 692mm (21.50" x 27.25") |
| 05.0109  | 546mm x 781mm (21.50" x 30.75") |
| 05.0110  | 629mm x 578mm (24.75" x 22.75") |



# Medical Device Packaging




## Coated Tyvek® Lid

### Features

Made from coated Tyvek®  
Sealing to thermoformed trays

Customized sizes and specifications supplied upon request



Authorized converter of DuPont™ Tyvek®  
Tyvek® is DuPont's trademark.



Lift Stickers



## Auto-Packing FFS - Form Fill Seal Packaging

### Features

Supplied with medical grade paper, Tyvek®, thermoforming film or other materials upon request

### Sterilization methods

Suitable for EtO and Irradiation



## Lift Stickers

### Features

Made from coated Tyvek®  
Stronger self-adhesives

Customized printing, sizes and specifications supplied upon request



# Medical Device Packaging



## Pre-made Pouches

### Features

- Made from transparent and milky film
- Chevron seal, thumb or tear notch,
- Tack seals and hang holes supplied upon request

### Sterilization methods

Suitable for Irradiation



## Breathable Strip Bags

### Features

- Made from transparent film with Tyvek® strip
- Transparent film on both sides for optimal content visibility and puncture-resistance
- Tear notch for easy open
- Outstanding and impervious medical sterilization packaging

### Sterilization methods

Suitable for EtO and Irradiation



## Aluminum Foil Bags

### Features

- Made from aluminum foil
- Chevron or corner seal
- Silver or white finishes
- Recommendation: Surgical blade, electrode

### Sterilization methods

Suitable for Irradiation



## Aluminum Foil / Tyvek® Header Bags

### Features

- Made from aluminum foil and Tyvek®
- Cleanliness and easy-to-open
- High sealing strength

### Sterilization methods

Suitable for EtO and Irradiation

Authorized converter of DuPont™ Tyvek®  
Tyvek® is DuPont's trademark.



# Medical Device Packaging

## Surgical Equipment Drape Packs

### Features

Made from transparent film  
Customized sizes and specifications supplied upon request

### Sterilization methods

Suitable for Irradiation



Robotic Arm Drape



Display Screen Drape



## Contact:

- ☎ **Office:** 886-2-2697-5499 **Fax:** 886-2-2697-5497  
8F-10, No. 99, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City 221, Taiwan
- Factory:** 886-2-2455-5498 **Fax:** 886-2-2455-5497  
1F, No. 33-7, Dahua 2nd RD., Qidu Dist., Keelung 206, Taiwan
- ✉ **E-mail:** [service@sterileright.com.tw](mailto:service@sterileright.com.tw)
- 🌐 **Website:** [www.sterileright.com.tw](http://www.sterileright.com.tw)