



SterileRight

Packaging Mfg. Inc.

www.sterileright.com.tw



Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.



Aseptic Transfer Packaging



EN868/ISO11607/11140/13485/ GMP



www.sterileright.com.tw

Steam Sterilization Bags

made from Tyvek® / film.

Produced in an ISO 5 (class 100) cleanroom



Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

- Made from **Tyvek®** and **HDPE** film
- Mainly made for **Steam** sterilization under controlled conditions
- The packed item can be used directly
- Intended usage for plastic or rubber bungs, plastic parts, or related products
- Steam** indicator ink printing is available
- Acceptable for **EtO** or **Gamma** sterilization as well.
- Customized sizes and specifications supplied upon request

Item No.	Sizes
11.7301	300mm x 500mm
11.7302	400mm x 500mm
11.7303	400mm x 600mm
11.7304	200mm x 300mm
11.7305	500mm x 650mm
11.7306	200mm x 350mm
11.7307	800mm x 800mm
11.7308	1000mm x 1000mm

01

EN868 ISO 11607 11140 13485






Easy Tear Steam Sterilization Bags

made from Tyvek® / film.

Produced in an ISO 5 (class 100) cleanroom



Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

- Made from **Tyvek®** and **HDPE** or **PP** film
- Mainly made for **Steam** sterilization under controlled conditions
- The packed item can be used directly
- Intended usage for plastic or rubber bungs, plastic parts, or related products
- Steam** indicator ink printing is available
- Customized sizes and specifications supplied upon request

Item No.	Sizes
15.0301	300mm x 500mm
15.0302	400mm x 500mm
15.0303	400mm x 600mm
15.0304	200mm x 300mm
15.0305	500mm x 650mm
15.0306	200mm x 350mm
15.0307	800mm x 800mm
15.0308	1000mm x 1000mm

02

EN868 ISO 11607 11140 13485



Steam Sterilization Rolls

made from Tyvek® / film.

Produced in an ISO 5 (class 100) cleanroom



Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

- Made from **Tyvek®** and **HDPE** film
- Mainly made for **Steam** sterilization under controlled conditions
- The packed item can be used directly
- Intended usage for plastic or rubber bungs, plastic parts, or related products
- Steam** indicator ink printing is available
- Acceptable for **EtO** or **Gamma** sterilization as well.
- Customized sizes and specifications supplied upon request

Item No.	Sizes
01.0505	200mm x 100M
01.0507	300mm x 100M
01.0509	400mm x 100M
01.0510	500mm x 100M
01.0511	600mm x 100M
01.0512	800mm x 100M
01.0513	1000mm x 100M

03



www.sterileright.com.tw

Easy Peel Steam Sterilization Pouches

made from Tyvek® / film.

Produced in an ISO 7 (class 10,000) cleanroom



Self-seal Steam Pouch



Peel Roll



Peel Pouch

Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

Made from **Tyvek®** /film

Mainly made for **Steam** sterilization under controlled conditions
The packed item can be used directly

Intended usage for goggles, garments, gloves, light plastic parts, or related products

Steam indicator ink printing is available

Customized sizes and specifications supplied upon request

Item No.

Sizes

04.0601	170mm x 265mm
04.0602	250mm x 375mm
04.0603	300mm x 460mm

Item No.

Sizes

01.0617	200mm x 100M
01.0618	300mm x 100M
01.0619	400mm x 100M

04

EN868 ISO 11607 11140 13485





www.sterileright.com.tw

Ultraclean Tyvek® Inserts (Sheets) made from Tyvek®

Produced in an ISO 5 (class 100) cleanroom



Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

Made from single Tyvek®

Intended primary usage as a clean protective layer for pre-fillable syringes

Customized sizes and specifications supplied upon request

05

EN868 ISO 11607 11140 13485





www.sterileright.com.tw

Clean LDPE Bags made from LDPE film.

Produced in an ISO 6 (class 1000) cleanroom



Pharmaceutical Use

Your reliable
OEM & ODM
partner

To market great and reliable
products worldwide.

Features

Made from **LDPE** film

Mainly made for **Gamma** sterilization

Clean LDPE Bag is the easy solution against contamination

Customized sizes and specifications supplied upon request

Item No.	Sizes
M901.02.0010-7	340mm x 510mm x 60um
M901.02.0010-1	400mm x 600mm x 60um
M901.02.0010-2	500mm x 800mm x 60um
M901.02.0010-3	500mm x 1000mm x 60um
M901.02.0010-4	500mm x 1200mm x 60um
M901.02.0010-5	600mm x 1000mm x 60um
M901.02.0010-6	600mm x 1200mm x 60um
M901.02.0010-8	340mm x 510mm x 100um
M901.02.0010-9	400mm x 600mm x 100um
M901.02.0010-10	500mm x 800mm x 100um
M901.02.0010-11	500mm x 1000mm x 100um
M901.02.0010-12	500mm x 1200mm x 100um
M901.02.0010-13	600mm x 1000mm x 100um
M901.02.0010-14	600mm x 1200mm x 100um

06

EN868 ISO 11607 11140 13485





Contact:

- 📍 **Office:** 886-2-2697-5499 **Fax:** 886-2-2697-5497
8F-10, No. 99, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City 221, Taiwan
- 🏭 **Factory:** 886-2-2455-5498 **Fax:** 886-2-2455-5497
1F, No. 33-7, Dahua 2nd RD., Qidu Dist., Keelung 206, Taiwan
- ✉️ **E-mail:** service@sterileright.com.tw
- 🌐 **Website:** www.sterileright.com.tw