



SterileRight

Packaging Mfg. Inc.



Medical Device **STERILIZATION PACKAGING**



EN868/ISO11607/11140/13485/ GMP www.sterileright.com.tw



Fairness Efficiency Innovation Development

- p04----- Reinforced Header Bags
- p05----- Peel Pouches
- p06----- FFS-Form Fill Seal Packaging
- p07----- Coated Tyvek® Lid
- p08----- Pre-made Pouches
Breathable Strip Bags
- p09----- Aluminum Foil Bags
Aluminum Foil / Tyvek® Header Bags
- p10----- Surgical Equipment Drape Packs



OEM & ODM



Reinforced Header Bags

Features

Made from un-coated or coated Tyvek® and high functional film
Ideal for heavy loading surgery and procedure kits
Suitable for bulk sterile kits, drapes and gowns
Higher seal strength & superior breathability
Impact strength, puncture and tear resistance

Customized sizes and specifications supplied upon request

Sterilization methods

Suitable for EtO and Irradiation



Patent in Taiwan

Medical Device Packaging



Peel Pouches

Features

Supplied with paper/paper, paper/film, Tyvek® /film and film/film
Chevron seal, thumb notch, with or without coating produced according to customer's need

Sterilization methods

Suitable for EtO



Sterilization Tyvek® Pouch



Customized printing, sizes and specifications supplied upon request

Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Item No. Sizes

05.0101	387mm x 425mm (15.25" x 16.75")
05.0102	413mm x 502mm (16.25" x 19.75")
05.0103	438mm x 578mm (17.25" x 22.75")
05.0104	464mm x 540mm (18.25" x 21.25")
05.0105	476mm x 413mm (18.75" x 16.25")
05.0106	514mm x 578mm (20.25" x 22.75")
05.0107	514mm x 629mm (20.25" x 24.75")
05.0108	546mm x 692mm (21.50" x 27.25")
05.0109	546mm x 781mm (21.50" x 30.75")
05.0110	629mm x 578mm (24.75" x 22.75")





OEM & ODM

Medical Device Packaging



Coated Tyvek® Lid


Features

Made from coated Tyvek®
Sealing to thermoformed trays

Customized sizes and specifications supplied upon request



Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.



Lift Stickers



FFS - Form Fill Seal Packaging

Features

Supplied with medical grade paper, Tyvek®, thermoforming film
or other materials upon request

Sterilization methods

Suitable for EtO and Irradiation



05

Lift Stickers

Features

Made from coated Tyvek®
Stronger self-adhesives

Customized printing, sizes and specifications supplied upon request



EN 868/ ISO 11607/11140/13485/ GMP

06

Medical Device Packaging



Pre-made Pouches

Features

Made from transparent and milky film
Chevron seal, thumb or tear notch,
Tack seals and hang holes supplied upon request

Sterilization methods

Suitable for Irradiation



Breathable Strip Bags

Features

Made from transparent film with Tyvek® strip
Transparent film on both sides for optimal content visibility
and puncture-resistance Tear notch for easy open
Outstanding and impervious medical sterilization packaging

Sterilization methods

Suitable for EtO and Irradiation



Aluminum Foil Bags

Features

Made from aluminum foil
Chevron or corner seal
Silver or white finishes
Recommendation: Surgical blade, electrode

Sterilization methods

Suitable for Irradiation



Aluminum Foil / Tyvek® Header Bags

Features

Made from aluminum foil and Tyvek®
Cleanliness and easy-to-open
High sealing strength

Sterilization methods

Suitable for EtO and Irradiation

Authorized converter of DuPont™ Tyvek®
Tyvek® is DuPont's trademark.

Medical Device Packaging

Surgical Equipment Drape Packs

Features

Made from transparent film
Customized sizes and specifications supplied upon request

Sterilization methods

Suitable for Irradiation



Robotic Arm Drape



Display Screen Drape





Contact:

- ☎ **Office:** 886-2-2697-5499 **Fax:** 886-2-2697-5497
8F-10, No. 99, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City 221, Taiwan
- Factory:** 886-2-2455-5498 **Fax:** 886-2-2455-5497
1F, No. 33-7, Dahua 2nd RD., Qidu Dist., Keelung 206, Taiwan
- ✉ **E-mail:** service@sterileright.com.tw
- 🌐 **Website:** www.sterileright.com.tw